

characteristics Labsit stool with castors

Article: 9127 integral foam

Dimensions in mm

| | | |
|-------|----|---------|
| Seat: | H: | 450-650 |
| | B: | 400 |
| | T: | 400 |

| | | | |
|--------------|--------|----------------|------|
| Packing data | Weight | kg | 7,4 |
| | Volume | m ³ | 0,19 |



Base equipment

Benefit

| | | |
|-----------------------|--|---|
| Laboratory compliance | Washable und resistant to disinfectant. GMP and cleanroom compatible. Can be used in laboratories with safety classification S1 to S3. Minimal joints and positive-locking (Hygienic-Design). | Suitable for laboratories and all hygiene-critical environments particularly health care, universities, biotechnology, chemistry and pharmaceuticals. |
| Upholstered Seat | Integral PU-foam The Integral PU-foam upholstery is extremely robust, durable, washable, resistant to disinfectant, PVC- free, resistant to blood and urine, fire protection according to EN 1021 Part 1 + 2 Colour: black. | For applications where aggressive substances are in use. The surface is hard-wearing, extremely robust, and easily cleaned. Resistant to pointed and sharp-edged equipment. |
| Seat Height | Seat height infinitely adjustable between 450 - 650 mm. | Large variability of seat height adjustment. |
| 5-Star Base | Sturdy polyamide base Optional: Aluminium base | Very sturdy, durable and stable. High quality. Attractive design. |
| Castors/Glides | Load dependant braked soft double castors for hard floors. | The chair will not roll away when unused which minimises the risk for accidents. |
| Standards | Suitable for cleanrooms category 3 in accordance with DIN EN ISO 14644-1 GMP-Conform Laboratory safety classification S1 to S3 Tested for harmful substances technical standard DIN 68 877-1 GS-Certified for highest levels of safety Blue Angel environmental certificate Qualitymanagement DIN EN ISO 9001 Environmental Management DIN EN ISO 14001 Verified Environmental Management EMAS Occupational Health & Safety OHSAS 18001 EcoVadis CSR Rating 2016 Gold. | Certified by independent testing institutes. Procurement security. |
| Guarantee | 10-year long-term-guarantee - 3 years of which full guarantee. | Quality assured and cost effective. |