

Pass-through refrigerator

PASS-LC 1400 GLASS

+2...+10°C



Standard equipment:

- ✓ RAL 9016 coated exterior
- ✓ Interior made of AISI 304 acid-resistant steel
- ✓ 70 mm insulation
- ✓ Environmentally friendly refrigerant
- ✓ 4 glass doors with replaceable seal
- ✓ Mechanical door lock
- ✓ Automatic door closing system`
- ✓ Electronic temperature controller with digital display with the possibility of changing the temperature setting every 0.1k
- ✓ Visual and sound alarm for high and low temperature
- ✓ Open door alarm
- ✓ Dirty condenser coil alarm
- ✓ Faulty temperature sensor alarm
- ✓ External alarm output potential-free contact
- ✓ Auto defrost and condensation evaporation
- ✓ Antifreeze system
- ✓ Shelf guides / horizontal supports
- ✓ 8 coated wire shelves
- ✓ Adjustable shelf position
- ✓ Sensor throughput
- ✓ 3 castors with brake and 2 castors without brake
- ✓ 2 additional adjustable feet
- ✓ Handles at the back of the unit for easy transport
- ✓ LED lighting

Optional equipment:

- ✓ Alarm backup battery
- ✓ Power failure alarm
- ✓ Temperature sensor
- ✓ Change of standard temperature range
- ✓ Additional control hole in the wall
- ✓ Glycol container
- ✓ Adjustable PVC feet
- ✓ 75 mm drawer with dividers
- ✓ 150 mm drawer
- ✓ Additional shelf guides / horizontal supports
- ✓ Additional shelves (coated wire, stainless steel wire)
- ✓ External housing made of acid-resistant AISI 304 steel

PASS-LC 1400 GLASS



Environmentally friendly refrigerant



Antifreeze system



Sensor throughput




2 castors with brake and 2 castors without brake



2 additional adjustable feet



LED lighting



Technical specifications:

	Width	Depth	Height
External dimensions [mm]	1480	870 (824)	1994
Internal dimensions [mm]	1340	680	1500
Cooling unit type	Internal		
Cooling type	Dynamic		
Gross capacity	1360 [dm ³]		
Net capacity	1260 [dm ³]		
Weight	205 [kg]		
Max. load	240 [kg]		
Number of shelves	8		
Number of guides / horizontal supports	16		
Max. shelf load	30 [kg]		
Temperature range	+2...+10 [°C]		
Ambient conditions	+25 [°C] op 60% HR		
Refrigerant	R290		
Power	230/50 [V/Hz]		
Number of external doors	2		
Daily energy consumption	8.0 [kWh/24h]		
Noise level	57 [db]		
Power rating	730 [W]		
Rated current	3.7 [A]		